


NEW JW-2 Water Pump



JW-2 Features

- Specially designed for endoscope examination and treatment.
- Operator can control water flow with foot switch.
- Large 1-liter water bottle is equipped for long time endoscope examination.
- Thanks to unique piping technology, operator can stop water flow very quickly in order to avoid water dropping.



Item	Specifications
Power supply	120V/230V
Power consumption rated	0.20A/0.12A
Fuse	T500mA/250Vx2
Applicable endoscopes	Fujinon endoscopes with a water-jet mechanism
Operation method	Turn on/off a foot switch
Flow rate adjustment	Unstepped continuously variable
Dimensions(WXHxD)	240X215X315mm(with standard accessories)
Mass	2.9kg

FUJIFILM

FUJIFILM Group

FUJINON Corporation

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA
331-9624, JAPAN
TEL : 81-48-668-2153 FAX : 81-48-668-1570
E-mail: sales@msv.fujinon.co.jp
URL <http://www.fujinon.co.jp/>

FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.
TEL : 1-800-490-0661 FAX : 1-973-633-8818
E-mail: med.info@fujinon.com
URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY
TEL : 49-2154-924-0 FAX : 49-2154-924-290
E-mail: Fujinon@fujinon.de
URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON
INDUSTRIAL PARK, SINGAPORE, 159552
TEL : 65-6276-4988 FAX : 65-6276-6911
E-mail: fujinon@fujinon.com.sg
URL <http://www.fujinon.com.sg/>





FUJINON AUSTRALIA PTY. LTD.




UNIT-18, 52 HOLKER STREET, SILVERWATER
N.S.W. 2128, AUSTRALIA
TEL : 61-2-9748-2744 FAX : 61-2-9748-2428
E-mail: sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD, 79 RI JING ROAD,
WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA
TEL : 86-21-5424-9262(Office) TEL : 86-21-5866-2151(Service)
FAX : 86-21-5866-2150
E-mail: fjng2001@hotmail.com

JW-2 Accessories

	Article no.	Quantity	Description
	68A11194990	1	Connection Adapter (for Scope)
	CT-1JW	2	Connecting Tube (Scope - Pump)
	FS-40JW	1	Foot Switch
	68A11188640	1	Maintenance Kit (Pump Tube+Grease)

	Article no.	Quantity	Description
	PU-1JW	1	Pump Cartridge
	WT-1000JW	1	Water Tank
	68A11915980	1	Drying Tube Adapter (for Maintenance)

FUJIFILM

FUJIFILM Group

FUJINON Corporation

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA
331-9624, JAPAN
TEL : 81-48-668-2153 FAX : 81-48-668-1570
E-mail: sales@msv.fujinon.co.jp
URL <http://www.fujinon.co.jp/>

FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.
TEL : 1-800-490-0661 FAX : 1-973-633-8818
E-mail: med.info@fujinon.com
URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILlich, GERMANY
TEL : 49-2154-924-0 FAX : 49-2154-924-290
E-mail: Fujinon@fujinon.de
URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON
INDUSTRIAL PARK, SINGAPORE, 159552
TEL : 65-6276-4988 FAX : 65-6276-6911
E-mail: fujinon@fujinon.com.sg
URL <http://www.fujinon.com.sg/>

FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER
N.S.W. 2128, AUSTRALIA
TEL : 61-2-9748-2744 FAX : 61-2-9748-2428
E-mail: sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD, 79 RI JING ROAD,
WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA
TEL : 86-21-5424-9262(Office) TEL : 86-21-5866-2151(Service)
FAX : 86-21-5866-2150
E-mail: fjng2001@hotmail.com